

Anti-Rbm15B Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ndt2

Item Number: CM0027862

Introduction

Rabbit polyclonal antibody against human RBM15B (Putative RNA-binding protein 15B, OTT3), suitable for WB and ELISA. Cross-reacts with human and mouse. Ideal for studying m6A methylation and RNA processing.

[Learn More](#)

Parameter	Value
Product Name	RBM15B Rabbit pAb
Remarks/Alias	OTT3; HsOTT3; HUMAGCGB; RBM15B
Species	Human
GeneID (Human)	29890
GeneID	29890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-720 of human RBM15B (NP_037418.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NDT2
Protein Weight	97kDa
Shipping	Ice bag