

Anti-Rbm17 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96I25

Item Number: CM0030293



Introduction

Rabbit polyclonal antibody raised against recombinant human RBM17 (SPF45) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects a 45 kDa nuclear splicing factor involved in mRNA splicing regulation.

[Learn More](#)

Parameter	Value
Product Name	RBM17 Rabbit pAb
Remarks/Alias	SPF45; RBM17
Species	Human
Gene ID (Human)	84991
Gene ID	84991
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RBM17 (NP_116294.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96I25
Protein Weight	45kDa
Shipping	Ice bag