

Anti-Rbm22 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nw64

Item Number: CM0016180



Introduction

High-quality rabbit polyclonal antibody targeting human RBM22 (Cwc2, ZC3H16, fSAP47). Validated for WB, ELISA; cross-reacts with human, mouse, rat. Ideal for spliceosome and pre-mRNA splicing research.

[Learn More](#)

Parameter	Value
Product Name	RBM22 Rabbit pAb
Remarks/Alias	Cwc2; ZC3H16; fSAP47; RBM22
Species	Human
Gene ID	55696 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human RBM22 (NP_060517.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NW64
Protein Weight	47kDa
Shipping	Ice bag