


Anti-Rbm25 Polyclonal Antibody For Wb, If/Icc, Elisa - P49756

Item Number: CM0018241



Introduction

Rabbit polyclonal antibody targeting RNA-binding protein 25 (RBM25), validated in WB, IF/ICC, ELISA; reacts with human and mouse. RBM25 regulates alternative splicing and apoptosis via BCL2L1 isoform modulation - P49756

[Learn More](#)

Parameter	Value
Product Name	RBM25 Rabbit pAb
Remarks/Alias	S164; NET52; RNPC7; Snu71; RED120; fSAP94; RBM25
Species	Human
GeneID (Human)	58517
GeneID	58517
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49756
Protein Weight	100kDa
Shipping	Ice bag