

Anti-Rbm26 Rabbit Polyclonal Antibody For Wb And Elisa - Q5T8P6

Item Number: CM0017760

Introduction

High-quality rabbit polyclonal antibody targeting RBM26, validated for WB and ELISA. Detects human RBM26 and cross-reacts with mouse and rat. Suitable for RNA-binding protein research and nuclear pA+ RNA turnover studies.

[Learn More](#)

Parameter	Value
Product Name	RBM26 Rabbit pAb
Remarks/Alias	ARRS2; SE70-2; ZC3H17; PRO1777; C13orf10; PPP1R132; RBM26
Species	Human
Gene ID (Human)	64062
Gene ID	64062
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5T8P6
Protein Weight	114kDa
Shipping	Ice bag