

Anti-Rbm26 Rabbit Polyclonal Antibody For Wb And Elisa - Q5T8P6

Item Number: CM0015375

Introduction

Anti-RBM26 rabbit polyclonal antibody applicable for Western blot and ELISA. Targets human RNA-binding protein 26, implicated in nuclear poly(A)+ RNA turnover. Ideal for post-transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	RBM26 Rabbit pAb
Remarks/Alias	ARRS2; SE70-2; ZC3H17; PRO1777; C13orf10; PPP1R132; RBM26
Species	Human
Gene ID (Human)	64062
Gene ID	64062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 65-140 of human RBM26 (NP_071401.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5T8P6
Protein Weight	114kDa
Shipping	Ice bag