

Anti-Rbm3 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P98179

Item Number: CM0002708



Introduction

Rabbit monoclonal antibody against human RBM3 (RNA-binding protein 3), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for cold stress and RNA-binding research.

[Learn More](#)

Parameter	Value
Product Name	RBM3 Rabbit mAb
Remarks/Alias	RNPL; IS1-RNPL; RBM3
Species	Human
GeneID (Human)	5935
GeneID	5935
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 58-157 of human RBM3 (P98179)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98179
Protein Weight	17kDa
Shipping	Ice bag