


Anti-Rbm39 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14498

Item Number: CM0019557



Introduction

High-quality rabbit polyclonal antibody against human RBM39 (CAPER α). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for splicing and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	RBM39 Rabbit pAb
Remarks/Alias	HCC1, CAPER, RNPC2, FSAP59, CAPERalpha, RBM39
Species	Human
GeneID (human)	9584
GeneID	9584
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-508 of human RBM39 (NP_001229528.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14498
Protein Weight	59kDa
Shipping	Ice bag