

Anti-Rbm46 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tby0

Item Number: CM0024129

Introduction

Rabbit polyclonal antibody targeting RBM46 (CT68), a key regulator of meiosis and fertility. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for reproductive biology and germ cell research.

[Learn More](#)

Parameter	Value
Product Name	RBM46 Rabbit pAb
Remarks/Alias	CT68; RBM46
Species	Human
GeneID (human)	166863
GeneID	166863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human RBM46 (NP_659416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBY0
Protein Weight	60kDa
Shipping	Ice bag