


Anti-Rbm5 Polyclonal Antibody For Wb, Elisa - P52756

Item Number: CM0014761



Introduction

Rabbit polyclonal antibody against human RBM5, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Suitable for apoptosis and alternative splicing research.

[Learn More](#)

Parameter	Value
Product Name	RBM5 Rabbit pAb
Remarks/Alias	G15, H37, RMB5, LUCA15, LUCA-15, RBM5
Species	Human
GeneID (human)	10181
GeneID	10181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RBM5 (NP_005769.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52756
Protein Weight	92kDa
Shipping	Ice bag