

Anti-Rbm6 Rabbit Polyclonal Antibody For Wb And Elisa - P78332

Item Number: CM0014299

Introduction

Rabbit polyclonal antibody against human RBM6, validated for Western blot and ELISA. Suitable for studying RNA-binding protein 6 function, alternative splicing, and lung cancer antigen research.

[Learn More](#)

Parameter	Value
Product Name	RBM6 Rabbit pAb
Remarks/Alias	3G2; g16; DEF3; DEF-3; HLC-11; NY-LU-12; RBM6
Species	Human
GeneID	10180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 530-770 of human RBM6 (NP_005768.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78332
Protein Weight	129kDa
Shipping	Ice bag