

Anti-Rbm7 (Ko Validated) Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - Q9Cqt2

Item Number: CM0010599

Introduction

KO validated rabbit monoclonal antibody against mouse Rbm7. Suitable for WB, IP, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Protein weight 30 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Rbm7 Rabbit mAb
Remarks/Alias	1200007M24Rik; 1500011D06Rik
Species	Mouse
GeneID (human)	10179
GeneID	67010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 87-265 of mouse Rbm7 (NP_659197.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9CQT2
Protein Weight	30kDa
Shipping	Ice bag