


Anti-Rbm7 Polyclonal Antibody For Wb, Elisa - Q9Y580

Item Number: CM0019655



Introduction

High-quality rabbit polyclonal antibody against human RNA-binding protein 7 (RBM7), validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for studying NEXT complex function and RNA metabolism.

[Learn More](#)

Parameter	Value
Product Name	RBM7 Rabbit pAb
Remarks/Alias	RBM7
Species	Human
GeneID (Human)	10179
GeneID	10179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of human RBM7 (NP_057174.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y580
Protein Weight	31kDa
Shipping	Ice bag