

Anti-RbmX/Hnrnp G Rabbit Polyclonal Antibody For Wb, Elisa - P38159

Item Number: CM0018372



Introduction

Rabbit polyclonal antibody recognizing human RBMX/hnRNP G, validated for WB and ELISA, cross-reacts with mouse. Targets a multifunctional RNA-binding protein involved in transcription, splicing, and TNFR1 signaling. Ideal for oncology and neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	RBMX/hnRNP G Rabbit pAb
Remarks/Alias	RNMX; HNRPG; HNRNPG; MRXS11; RBMXP1; RBMXRT; hnRNP-G; RBMX/hnRNP G
Species	Human
GeneID (Human)	27316
GeneID	27316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human RBMX/hnRNP G (NP_002130.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P38159
Protein Weight	42kDa
Shipping	Ice bag