

Anti-Rbmy1F Polyclonal Antibody For Wb And Elisa - Q15415

Item Number: CM0018849

Introduction

High-quality rabbit polyclonal antibody targeting human RBMY1F (UniProt Q15415). Suitable for Western blot and ELISA applications, with rat cross-reactivity. Essential for studying spermatogenesis and pre-mRNA splicing in testis.

[Learn More](#)

Parameter	Value
Product Name	RBMY1F Rabbit pAb
Remarks/Alias	YRRM2; RBMY1F
Species	Human
GeneID (Human)	159163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human RBMY1F (NP_689798.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q15415
Protein Weight	56kDa
Shipping	Ice bag