


Anti-Rbp1 Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P09455

Item Number: CM0014546



Introduction

High-quality rabbit polyclonal antibody targeting human RBP1 (cellular retinol-binding protein I). Validated for WB, IF/ICC, IF-P, ELISA; cross-reacts with mouse and rat. Ideal for vitamin A transport and lipid droplet biology studies.

[Learn More](#)

Parameter	Value
Product Name	RBP1 Rabbit pAb
Remarks/Alias	CRBP; RBPC; CRBP1; CRBPI; hCRBP1; CRABP-I; RBP1
Species	Human
GeneID (human)	5947
GeneID	5947
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RBP1 (NP_001352869.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09455
Protein Weight	16kDa
Shipping	Ice bag