

Anti-Rbp3 Rabbit Polyclonal Antibody For Wb And Elisa - P10745

Item Number: CM0025362

Introduction

Rabbit polyclonal antibody against human RBP3, validated for Western blot and ELISA applications. Detects human, mouse, and rat RBP3 (135 kDa), a key interphotoreceptor retinoid-binding protein in the visual cycle.

[Learn More](#)

Parameter	Value
Product Name	RBP3 Rabbit pAb
Remarks/Alias	IRBP; RBPI; RP66; D10S64; D10S65; D10S66; RBP3
Species	Human
GeneID	5949
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-320 of human RBP3 (NP_002891.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10745
Protein Weight	135kDa
Shipping	Ice bag