


Anti-Rbp4 Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P02753

Item Number: CM0010936



Introduction

RBP4 Rabbit Monoclonal Antibody (alias RDCCAS; MCOPCB10) validated for WB, IF/ICC, IHC-P, and ELISA. Reacts with human, mouse, and rat retinol-binding protein 4 (plasma retinol transport). Molecular weight ~21 kDa.

[Learn More](#)

Parameter	Value
Product Name	RBP4 Rabbit mAb
Remarks/Alias	RDCCAS; MCOPCB10
Species	Human
GeneID (Human)	5950
GeneID	5950
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human RBP4 (NP_006735.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02753
Protein Weight	21kDa
Shipping	Ice bag