


## Anti-Rbp4 Polyclonal Antibody For Wb, Elisa - P02753

Item Number: CM0015617



### Introduction

Rabbit polyclonal antibody targeting human RBP4 (retinol-binding protein 4), validated for WB and ELISA applications. Detects endogenous levels across human, mouse, and rat species. Ideal for retinol transport and vitamin A metabolism research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | RBP4 Rabbit pAb                                                                                                  |
| Remarks/Alias    | RDCCAS; MCOPCB10; RBP4                                                                                           |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 5950                                                                                                             |
| Gene ID          | 5950                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human RBP4 (NP_006735.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P02753                                                                                                           |
| Protein Weight   | 21kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |