


Anti-Rbp4 Rabbit Polyclonal Antibody For Wb, Elisa - P02753

Item Number: CM0021829



Introduction

High-quality rabbit polyclonal antibody targeting human RBP4, suitable for Western blot and ELISA. Cross-reacts with mouse. Ideal for investigations into retinol transport, vitamin A homeostasis, and vision-related pathways.

[Learn More](#)

Parameter	Value
Product Name	RBP4 Rabbit pAb
Remarks/Alias	RDCCAS; MCOPCB10; RBP4
Species	Human
GeneID (Human)	5950
GeneID	5950
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human RBP4 (NP_006735.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02753
Protein Weight	21kDa
Shipping	Ice bag