


Anti-Rbp7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96R05

Item Number: CM0024327



Introduction

Rabbit polyclonal anti-RBP7 antibody for WB, IF/ICC, ELISA. Recognizes human, mouse, rat RBP7 (CRBP4/CRABP4). Key intracellular retinol transport protein. Ideal for retinoid metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RBP7 Rabbit pAb
Remarks/Alias	CRBP4; CRABP4; CRBP4V; RBP7
Species	Human
GeneID (Human)	116362
GeneID	116362
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human RBP7 (NP_443192.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96R05
Protein Weight	16kDa
Shipping	Ice bag