

Anti-Rbpjk Monoclonal Antibody For Wb, Ihc-P, Elisa - Q06330

Item Number: CM0005467



Introduction

Rabbit monoclonal antibody against RBPJK (CBF-1, SUH) validated for WB, IHC-P, and ELISA. Reacts with human and mouse RBPJK (56 kDa). Essential tool for studying Notch signaling and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	RBPJK Rabbit mAb
Remarks/Alias	SUH; csl; AOS3; CBF1; KBF2; CBF-1; RBP-J; RBPJK; IGKJRB; RBP-JK; RBPSUH; IGKJRB1; RBP-J kappa
Species	Human
GeneID (human)	3516
GeneID	3516
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RBPJK (Q06330).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06330
Protein Weight	56kDa
Shipping	Ice bag