

Anti-Rbpjk Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q06330

Item Number: CM0024989



Introduction

High-quality rabbit polyclonal antibody against human RBPJK (RBP-J, CBF1). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for Notch signaling and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	RBPJK Rabbit pAb
Remarks/Alias	SUH; csl; AOS3; CBF1; KBF2; CBF-1; RBP-J; RBPJK; IGKJRB; RBP-JK; RBPSUH; IGKJRB1; RBP-J kappa
Species	Human
GeneID (Human)	3516
GeneID	3516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human RBPJKB (NP_976028.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06330
Protein Weight	56 kDa
Shipping	Ice bag