

Anti-Rbpms Rabbit Polyclonal Antibody For Wb And Elisa - Q93062

Item Number: CM0016285



Introduction

Rabbit polyclonal antibody targeting RBPMS, an RNA-binding protein involved in alternative splicing regulation. Validated for WB and ELISA in human, mouse, and rat. Ideal for vascular biology and splicing research.

[Learn More](#)

Parameter	Value
Product Name	RBPMS Rabbit pAb
Remarks/Alias	HERMES; RBPMS
Species	Human
Gene ID	11030
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-219 of human RBPMS (NP_001008712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93062
Protein Weight	22kDa
Shipping	Ice bag