

Anti-Rc3H1/Roquin-1 Polyclonal Antibody For Wb & Elisa - Q5Tc82

Item Number: CM0019225



Introduction

High-quality rabbit polyclonal antibody targeting human RC3H1 (Roquin-1), validated for WB and ELISA. Suitable for studying post-transcriptional gene regulation and immune responses. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	RC3H1 Rabbit pAb
Remarks/Alias	FHL6; IMDSHY; RNF198; ROQUIN; RC3H1
Species	Human
Gene ID (Human)	149041
Gene ID	149041
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 834-1133 of human RC3H1 (NP_742068.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5TC82
Protein Weight	126kDa
Shipping	Ice bag