

Anti-Rcc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18754

Item Number: CM0027674



Introduction

Rabbit polyclonal antibody targeting RCC1, validated for WB, IHC-P, IF/ICC, and ELISA across human, mouse, and rat samples. Essential for studies on chromosome condensation and Ran-mediated nuclear processes.

[Learn More](#)

Parameter	Value
Product Name	RCC1 Rabbit pAb
Remarks/Alias	CHC1; RCC1-I; RCC1
Species	Human
Gene ID (Human)	1104
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RCC1 (NP_001260.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18754
Protein Weight	45kDa
Shipping	Ice bag