

Anti-Rcc1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P18754

Item Number: CM0002948



Introduction

RCC1 Rabbit mAb is a monoclonal antibody raised against human RCC1, validated for WB, IHC-P, and ELISA, with cross-reactivity in human, mouse, and rat. Ideal for studying chromatin condensation and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	RCC1 Rabbit mAb
Remarks/Alias	CHC1; RCC1-I; RCC1
Species	Human
GeneID (human)	1104
GeneID	1104
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RCC1 (P18754).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18754
Protein Weight	45kDa
Shipping	Ice bag