

Anti-Rcc2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P258

Item Number: CM0016250



Introduction

Rabbit polyclonal antibody targeting human RCC2 (TD-60) protein, validated for Western blotting and ELISA. Ideal for investigating cell cycle regulation, mitosis, and small GTPase signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	RCC2 Rabbit pAb
Remarks/Alias	TD-60; RCC2
Species	Human
Gene ID (Human)	55920
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 323-522 of human RCC2 (NP_001129676.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9P258
Protein Weight	56kDa
Shipping	Ice bag