


Anti-Rcl1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2P8

Item Number: CM0025636



Introduction

High-quality rabbit polyclonal antibody targeting human RCL1 (RNA 3'-terminal phosphate cyclase-like protein). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for ribosomal biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	RCL1 Rabbit pAb
Remarks/Alias	RNAC; RPCL1; RCL1
Species	Human
Gene ID	10171
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RCL1 (NP_005763.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2P8
Protein Weight	41kDa
Shipping	Ice bag