

Anti-Rcn1 Rabbit Polyclonal Antibody For Wb/If/Icc/Elisa - Q15293

Item Number: CM0018863

Introduction

High-quality rabbit polyclonal antibody against human RCN1 (Reticulocalbin-1), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for ER calcium signaling research.

[Learn More](#)

Parameter	Value
Product Name	RCN1 Rabbit pAb
Remarks/Alias	RCN; RCAL; PIG20; HEL-S-84; RCN1
Species	Human
GeneID (human)	5954
GeneID	5954
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-331 of human RCN1 (NP_002892.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15293
Protein Weight	39kDa
Shipping	Ice bag