

Anti-Rcn2 (Reticulocalbin-2) Polyclonal Antibody For Wb And Elisa - Q14257

Item Number: CM0018939

Introduction

Rabbit polyclonal antibody against human RCN2 (Reticulocalbin-2), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for endoplasmic reticulum calcium-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RCN2 Rabbit pAb
Remarks/Alias	E6BP; ERC55; ERC-55; TCBP49; N2
Species	Human
Gene ID	5955
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-180 of human RCN2 (NP_002893.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14257
Protein Weight	37kDa
Shipping	Ice bag