


Anti-Rcvrn Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P35243

Item Number: CM0029325



Introduction

High-quality rabbit polyclonal antibody targeting human recoverin (RCVRN). Validated for WB, IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for visual transduction research and cancer-associated retinopathy studies.

[Learn More](#)

Parameter	Value
Product Name	RCVRN Rabbit pAb
Remarks/Alias	RCV1; RCVRN
Species	Human
Gene ID (Human)	5957
Gene ID	5957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human RCVRN (NP_002894.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35243
Protein Weight	23 kDa
Shipping	Ice bag