

Anti-Rdh11 Polyclonal Antibody For Wb And Elisa - Q8Tc12

Item Number: CM0021995

Introduction

Anti-RDH11 polyclonal antibody suitable for WB and ELISA, detecting human and rat retinol dehydrogenase 11. Targets an enzyme critical in retinoid metabolism with a predicted molecular weight of 35 kDa. Ideal for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	RDH11 Rabbit pAb
Remarks/Alias	MDT1; CGI82; PSDR1; RALR1; SCALD; ARSDR1; HCBP12; RDJCSS; SDR7C1; RDH11
Species	Human
Gene ID (Human)	51109
Gene ID	51109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-145 of human RDH11 (NP_057110.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8TC12
Protein Weight	35kDa
Shipping	Ice bag