

Anti-Rdm1 (Rad52B) Polyclonal Antibody For Wb, Elisa - Q8Ng50

Item Number: CM0030558

Introduction

RDM1 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with mouse. Targets human RAD52 motif-containing protein 1, a DNA/RNA-binding protein linked to cisplatin resistance.

[Learn More](#)

Parameter	Value
Product Name	RDM1 Rabbit pAb
Remarks/Alias	RAD52B; RDM1
Species	Human
GenelD (Human)	201299
GenelD	201299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-284 of human RDM1 (NP_663629.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8NG50
Protein Weight	32kDa
Shipping	Ice bag