

Anti-Reck Polyclonal Antibody For Wb, Elisa - 095980

Item Number: CM0028369



Introduction

Rabbit polyclonal antibody targeting human RECK (reversion-inducing cysteine-rich protein with Kazal motifs). Suitable for Western blot and ELISA, cross-reacts with mouse. Ideal for Wnt signaling, angiogenesis, and tumor biology studies.

[Learn More](#)

Parameter	Value
Product Name	RECK Rabbit pAb
Remarks/Alias	ST15; RECK
Species	Human
GeneID (Human)	8434
GeneID	8434
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 110-210 of human RECK (NP_066934.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95980
Protein Weight	106kDa
Shipping	Ice bag