

Anti-Recql4 Rabbit Polyclonal Antibody For Wb And Elisa - O94761

Item Number: CM0026678



Introduction

Rabbit polyclonal antibody targeting human RECQL4, validated for Western blot and ELISA. Detects RECQL4 in human and mouse samples. Suitable for studying DNA helicase function and genomic stability.

[Learn More](#)

Parameter	Value
Product Name	RECQL4 Rabbit pAb
Remarks/Alias	RECQ4; RECQL4
Species	Human
Gene ID (Human)	9401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 231-330 of human RECQL4 (NP_004251.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94761
Protein Weight	133kDa
Shipping	Ice bag