

Anti-Reep1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9H902

Item Number: CM0026892

Introduction

Rabbit polyclonal antibody targeting human REEP1 (Spastic paraplegia 31 protein), validated for WB and ELISA. Detects human, mouse, and rat REEP1 (~22 kDa). Suitable for ER network and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	REEP1 Rabbit pAb
Remarks/Alias	DSMA6; HMN5B; SPG31; Yip2a; C2orf23; REEP1
Species	Human
Gene ID (Human)	65055
Gene ID	65055
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H902
Protein Weight	22kDa
Shipping	Ice bag