

## Anti-Reep2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Brk0

Item Number: CM0021825

### Introduction

REEP2 rabbit polyclonal antibody for WB, IF-P, ELISA. Detects human REEP2, cross-reactive with mouse and rat. Targets ER network formation protein implicated in hereditary spastic paraplegia. Ideal for neuroscience and cell biology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | REEP2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SPG72; Yip2d; C5orf19; SGC32445; REEP2                                                                             |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 51308                                                                                                              |
| Gene ID          | 51308                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 93-252 of human REEP2 (NP_057690.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IF-P, ELISA                                                                                                    |
| Cross Reactivity | Mouse, Rat                                                                                                         |
| SWISS            | Q9BRK0                                                                                                             |
| Protein Weight   | 28kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |