

Anti-Reep5 Polyclonal Antibody For Wb, Elisa - Q00765

Item Number: CM0016795

Introduction

Rabbit polyclonal antibody against REEP5, validated for Western blot (WB) and ELISA; detects human, mouse, rat REEP5. Ideal for studying cardiac sarcoplasmic reticulum function. Immunogen corresponds to residues 120-189.

[Learn More](#)

Parameter	Value
Product Name	REEP5 Rabbit pAb
Remarks/Alias	DP1; TB2; YOP1; POB16; Yip2e; D5S346; C5orf18; REEP5
Species	Human
GeneID (Human)	7905
GeneID	7905
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-189 of human REEP5 (NP_005660.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00765
Protein Weight	21kDa
Shipping	Ice bag