


Anti-Reg1A Polyclonal Antibody For Wb And Elisa - P05451

Item Number: CM0029465



Introduction

High-quality rabbit polyclonal antibody against human REG1A (Lithostathine-1-alpha), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for pancreatic and brain regeneration research.

[Learn More](#)

Parameter	Value
Product Name	REG1A Rabbit pAb
Remarks/Alias	P19, PSP, PTP, REG, ICRF, PSPS, PSPS1, REG1A
Species	Human
GeneID (Human)	5967
GeneID	5967
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human REG1A (NP_002900.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05451
Protein Weight	19kDa
Shipping	Ice bag