


Anti-Reg3A Rabbit Polyclonal Antibody For Wb, Elisa - Q06141

Item Number: CM0019808



Introduction

Validate REG3A expression with this rabbit polyclonal antibody, suitable for WB and ELISA applications. Detects human, mouse, and rat REG3A, a bactericidal C-type lectin involved in inflammation and tissue regeneration.

[Learn More](#)

Parameter	Value
Product Name	REG3A Rabbit pAb
Remarks/Alias	HIP, PAP, PAP1, REG3, INGAP, PAP-H, PBCGF, HIP/PAP, REG-III, REG3A
Species	Human
Gene ID	5068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-175 of human REG3A (NP_002571.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06141
Protein Weight	19kDa
Shipping	Ice bag