


Anti-Reg3A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q06141

Item Number: CM0014712



Introduction

High-quality Anti-REG3A Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Cross-reacts with Human, Mouse, Rat. Targets REG3A (Q06141) for antimicrobial and inflammatory response studies.

[Learn More](#)

Parameter	Value
Product Name	REG3A Rabbit pAb
Alias	HIP, PAP, PAP1, REG3, INGAP, PAP-H, PBCGF, HIP/PAP, REG-III, REG3A
Species	Human
Gene ID	5068
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-175 of human REG3A (NP_002571.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06141
Protein Weight	19kDa
Shipping	Ice bag