

Anti-Reg4 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Byz8

Item Number: CM0013344



Introduction

Rabbit polyclonal antibody against human REG4, suitable for WB and ELISA. Cross-reacts with mouse and rat. Targets regenerating islet-derived protein 4, a lectin involved in gastrointestinal inflammation and cancer. Ideal for colorectal cancer research.

[Learn More](#)

Parameter	Value
Product Name	REG4 Rabbit pAb
Remarks/Alias	GISP; RELP; REG-IV; REG4
Species	Human
GeneID	83998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-158 of human REG4 (NP_114433.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BYZ8
Protein Weight	18kDa
Shipping	Ice bag