

Anti-Reln Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P78509

Item Number: CM0024469



Introduction

RELN Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat RELN. Recombinant immunogen (aa 50-190). For neuroscience and developmental biology research.

[Learn More](#)

Parameter	Specification
Product Name	RELN Rabbit pAb
Remarks/Alias	RL; ETL7; LIS2; PRO1598; RELN
Species	Human
Gene ID (Human)	5649
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-190 of human RELN (NP_005036.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78509
Protein Weight	388kDa
Shipping	Ice bag