

Anti-Renin Monoclonal Antibody For Wb, If/Icc, Elisa - P00797

Item Number: CM0004992



Introduction

Rabbit monoclonal antibody against human Renin (P00797). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying the renin-angiotensin system.

[Learn More](#)

Parameter	Value
Product Name	Renin Rabbit mAb
Remarks/Alias	RTD; HNFJ2; ADTKD4; Renin
Species	Human
Gene ID (Human)	5972
Gene ID	5972
Immunogen	Recombinant protein corresponding to a sequence within amino acids 250-406AA of human Renin (NP_000528.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00797
Protein Weight	45kDa
Shipping	Ice bag