

Anti-Renin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P00797

Item Number: CM0021697



Introduction

Renin Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, and rat renin. A research tool for renin-angiotensin system studies, hypertension, and kidney function analysis.

[Learn More](#)

Parameter	Value
Product Name	Renin Rabbit pAb
Remarks/Alias	RTD; HNFJ2; ADTKD4; Renin
Species	Human
Gene ID (Human)	5972
Gene ID	5972
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-406 of human Renin (NP_000528.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00797
Protein Weight	45kDa
Shipping	Ice bag