

Anti-Repin1 Polyclonal Antibody For Wb, Elisa - Q9Bwe0

Item Number: CM0023589

Introduction

REPIN1 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human and mouse REPIN1, a 64 kDa DNA-binding protein involved in lipid and iron metabolism regulation. Ideal for cell biology and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	REPIN1 Rabbit pAb
Remarks/Alias	AP4; RIP60; ZNF464; Zfp464; REPIN1
Species	Human
Gene ID	29803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 376-567 of human REPIN1 (NP_037532.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BWE0
Protein Weight	64kDa
Shipping	Ice bag