

## Anti-Repin1 Rabbit Monoclonal Antibody For Wb - Q9Bwe0

Item Number: CM0005059

### Introduction

Rabbit monoclonal REPIN1 antibody validated for WB, IHC-P, ELISA. Recognizes the 60 kDa DNA-binding protein involved in fatty acid import and lipid metabolism, ideal for metabolic research and IVD assay development.

[Learn More](#)

| Parameter        | Specification                                                                                                       |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | REPIN1 Rabbit mAb                                                                                                   |
| Remarks/Alias    | AP4; RIP60; ZNF464; Zfp464                                                                                          |
| Species          | Human                                                                                                               |
| GeneID (human)   | 29803                                                                                                               |
| GeneID           | 29803                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 376-567 of human REPIN1 (NP_037532.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | Q9BWE0                                                                                                              |
| Protein Weight   | 64kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |