

Anti-Reptin/Ruvbl2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y230

Item Number: CM0011727



Introduction

High-quality rabbit monoclonal antibody against human Reptin/RUVBL2 (RUVBL2, TIP48, ECP51). Validated for WB and ELISA, cross-reacts with mouse. Ideal for chromatin remodeling, transcription, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Reptin/RUVBL2 Rabbit mAb
Remarks/Alias	RVB2; TIH2; ECP51; TIP48; CGI-46; ECP-51; INO80J; REPTIN; TIP49B; TAP54-beta; Reptin/RUVBL2
Species	Human
GeneID (Human)	10856
GeneID	10856
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Reptin/RUVBL2 (Q9Y230)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y230
Protein Weight	51kDa
Shipping	Ice bag