

Anti-Rere (Atrophin-1-Like Protein) Polyclonal Antibody For Wb And Elisa - Q9P2R6

Item Number: CM0013180



Introduction

Rabbit polyclonal antibody targeting human RERE (Atrophin-1-like protein), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying transcriptional repression and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	RERE Rabbit pAb
Remarks / Alias	ARG, ARP, DNB1, ATN1L, NEDBEH, RERE
Species	Human
GeneID	473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human RERE (NP_036234.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2R6
Protein Weight	172kDa
Shipping	Ice bag